

Get a firm grip on sensitivity.

Ultraderm® gloves are ideal for those who require a glove material with reduced concentrations of chemical accelerators that are known to be potential sensitizers. These gloves have reduced amounts of natural rubber latex proteins.

Ultraderm gloves with enhanced comfort are constructed with a three-plane mold design and tapered fingers. The anti-slip finish offers excellent grip properties when handling instruments and the tapered cuff helps keep the glove firmly in place. Packaged in the smaller, easy-to-store Allegiance pack.



Ordering Information

Catalog Number	Size
2D7260	5.5
2D7261	6
2D7262	6.5
2D7263	7
2D7264	7.5
2D7265	8
2D7266	8.5
2D7267	9

All products packaged
40 pairs per box / 5 boxes per case

Professional Services at Allegiance

Our Professional Services Department offers ongoing service and support to ensure customer satisfaction. Our toll-free number provides a direct source for inquiries regarding clinical application and product usage. The Professional Services Group is staffed by registered nurses who offer a full range of clinical support services and product information.

To contact the Professional Services Department, call 1-800-327-7503.

- Reduced concentrations of chemical accelerators that are known to be potential sensitizers.
- These latex gloves contain 50 micrograms or less of total water-extractable protein per gram. *Caution: Safe use of this glove by or on latex sensitized individuals has not been established.*
- Exceed ASTM requirements for freedom from holes.
- Breakthrough time of glutaraldehyde 3.56 hours using a modification of ASTM F1407-96.
- Full 12" length for added protection.
- 8.7 mil average finger thickness.
- Excellent grip properties for handling instruments.
- Patented Triflex® mold design with tapered fingers.
- Tapered cuff prevents rolldown and helps glove stay in place.
- Packaged in smaller, easy to open Allegiance pack.

To order, call Allegiance Customer Service at 1-800-964-5227.



a Cardinal Health company

**Allegiance Healthcare Corporation
Gloves**

1500 Waukegan Road
McGaw Park, Illinois 60085

www.cardinal.com/allegiance

©Copyright 2000, Allegiance Healthcare Corporation. All rights reserved.
Lit. No. GLV00620 (0900/5M)

♻️ Printed on recycled paper with 10% post-consumer fiber.

ULTRADERM

**Reduced
proteins**

**Fewer potential
sensitizers**

**Excellent
grip**

**Enhanced
comfort**

**Outstanding
protection**

Get a Firm Grip on Sensitivity.